


Peel & Stick Label Options

Perfect For Barcoding Any Non-Textile





2" X 1" custom order



1.11" x .67" custom order



 $1^{1}/_{2}$ " $\times 3/_{4}$ " custom order



 $2^{1/2}$ " x 1" custom order



(Labels shown actual size)

Description: 8-digit, code interleaved 2 of 5 stocked. Different codes & sizes can be printed to

customer specifications.

Stock Size: 1.5" x .5" (Generic labels in stock for immediate delivery.)

Material: White Polyester (Over Laminated With Clear Polyester)

Adhesive: Aggressive Adhesive

Labels are supplied in rolls of 1,000. Simply peel & stick! Labels stick to most surfaces.

Shoe Recommendations: For easy scanning we recommend placement in the arch or side panel area of the shoe. Clean the area with alcohol prior to label placement for best adhesion. This should also be done on new merchandise to remove any residue that may have been deposited during manufacture. Allow alcohol to dry completely prior to label placement. After placement, labels should sit 24 hours prior to use to allow proper curing of the adhesive. Higher initial bonds can be achieved through increased rubdown pressure. Depending on how you maintain your inventory, you may barcode one shoe or both with either a unique or identical serial number.

For ordering labels, prices, or more information contact Amber@ipc-tags.com Customer Service: 817-714-8400 • Order Fax: (817) 887-1725 • Corporate (972) 822-0934